



NOVOSYN[®]


ABSORBABLE SUTURES

NOVOSYN[®]

ABSORBABLE SUTURES

Description

Novosyn[®] is a mid-term absorbable braided and coated synthetic suture made of polyglactin 910. Novosyn[®] degrades by hydrolysis and provides a predictable and reliable absorption.



Novosyn[®]




Product overview

Polyglactin 910, braided, coated



Color	Violet or undyed
Origin	Synthetic
Sizes	USP 8/0 (metric 0.4) to USP 2 (metric 5)
Type of absorption	By hydrolysis
Mass absorption	From 56 to 70 days
Sterilisation	Ethylene oxide








	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
Pre-cut sutures										
5 threads of 70 cm	5 x 70 violet						B0058449	B0058450		
10 threads of 60 cm	10 x 60 violet					B0058567	B0058568			
6 threads of 45 cm	6 x 45 violet					B0058647	B0058649	B0058650		
12 threads of 45 cm	12 x 45 violet					B0058667	B0058669	B0058670		
1 thread of 150 cm	1 x 150 violet		B0058704	B0058705	B0058706	B0058707	B0058709	B0058708		
5 threads of 70 cm	5 x 70 undyed								B0059450	
6 threads of 45 cm	6 x 45 undyed					B0059647	B0059649			
12 threads of 45 cm	12 x 45 undyed					B0059667	B0059669	B0059670		







Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle										
 HR 65	70 70 violet			B0068071		B0068072				
 HR 76	70 violet		B0068073	B0068074	B0068075	B0068076	B0068077			
 HR 48	150 violet					B0068119				

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

NOVOSYN[®]







ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle with short cutting point										
 HRC 37 Easyslide	4 x 70 violet					B0068124				
1/2 circle round bodied needle										
 HR 37s Easyslide	150 violet						B0068125			
1/2 circle round bodied needle with short cutting point										
 HRC 37 Easyslide	90 violet						B0068143			
1/2 circle round bodied needle										
 HR 70 Easyslide	70 violet					B0068183				
1/2 circle round bodied needle with blunt point										
 HRN 50s Easyslide	90 violet						B0068229			
1/2 circle reverse cutting needle										
 HS 37s + HR 37s Easyslide	150 violet			B0068278						
 HS 48 Easyslide	150 violet					B0068303				

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 HRC 37ss	90 violet				B0068344				
1/2 circle round bodied needle with trocar point									
 HRT 43	150 violet				B0068346	B0068347			
1/2 circle round bodied needle									
 HR 40s	150 violet					B0068350			
1/2 circle round bodied needle with short cutting point									
 HRC 37ss	90 violet					B0068354			
1/2 circle round bodied needle									
 HR 40s	150 violet				B0068355				
1/2 circle round bodied needle with trocar point									
 HRT 43s	150 violet				B0068356	B0068357			

NOVOSYN®








ABSORBABLE SUTURES








Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 <p>HRC 37ss</p> <p>Easyslide</p>	90 violet					B0068358			
1/2 circle round bodied needle with trocar point									
 <p>HRT 48</p> <p>Easyslide</p>	150 violet				B0068369	B0068370			
1/2 circle round bodied needle									
 <p>HR 60</p> <p>Easyslide</p>	150 violet					B0068380			
1/2 circle round bodied needle with short cutting point									
 <p>HRC 43ss</p> <p>Easyslide</p>	90 90 violet					B0068399	B0068404		
3/8 circle round bodied needle									
 <p>DR 48</p> <p>Easyslide</p>	90 violet					B0068437			
1/2 circle reverse cutting needle									
 <p>HS 48</p> <p>Easyslide</p>	90 violet					B0068463			

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with trocar point									
	90 violet							B0068473	
1/2 circle reverse cutting needle									
	90 violet							B0068488	
	90 violet							B0068498	
1/2 circle round bodied needle with trocar point									
	90 150 violet						B0068520	B0068503	
	90 violet							B0068535	
1/2 circle round bodied needle									
	120 90 violet						B0068555	B0068558	
	90 violet							B0068563	

NOVOSYN[®]







ABSORBABLE SUTURES







Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle with short cutting point										
 HRC 48 Easyslide	90 violet						B0068575			
1/2 circle round bodied needle										
 HR 76 Easyslide	90 violet					B0068579				
1/2 circle reverse cutting needle										
 HS 40s Easyslide	90 violet						B0068588			
1/2 circle round bodied needle										
 HR 37s Easyslide	90 violet						B0068598			
1/2 circle reverse cutting needle										
 HS 37s Easyslide	90 violet						B0068599			
1/2 circle round bodied needle with short cutting point										
 HRC 43s Easyslide	150 violet				B0068601	B0068602				
3/8 circle round bodied needle with blunt point										
 DRN 64 Easyslide	90 violet				B0068607	B0068605				

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle reverse cutting needle									
 HS 37s Easyslide	120 violet				B0068611				
3/8 circle round bodied needle with blunt point									
 DRN 64 Easyslide	90 violet						B0068614		
1/2 circle reverse cutting needle									
 HS 37s Easyslide	120 violet					B0068632			
1/2 circle round bodied needle									
 HR 65 Easyslide	90 violet				B0068639	B0068640			
1/2 circle round bodied needle with blunt point									
 HRN 45 Easyslide	90 violet						B0068648		
1/2 circle round bodied needle with trocar point									
 HRT 37 Easyslide	90 violet						B0068654		
1/2 circle round bodied needle									
 HR 37ss Easyslide	90 violet						B0068655		

NOVOSYN[®]








ABSORBABLE SUTURES








Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle										
 HR 40 Easyslide	90 violet						B0068661			
1/2 circle round bodied needle with blunt point										
 HRN 65 Easyslide	90 violet						B0068674	B0068675		
1/2 circle round bodied needle										
 HR 48 Easyslide	120 violet								B0068779	
 HR 70 Easyslide	70 violet				B0068803	B0068804				
1/2 circle round bodied needle with short cutting point										
 HRC 48 Easyslide	150 violet						B0068837			
1/2 circle round bodied needle										
 HR 43 Easyslide	90 violet								B0068848	

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 <p>HRC 43</p> <p>150 violet</p> <p>B0068872 B0068873</p>									
 <p>HRC 37s</p> <p>90 violet</p> <p>B0068897</p>									
1/2 circle round bodied needle with trocar point									
 <p>HRT 48</p> <p>150 violet</p> <p>B0068918</p>									
1/2 circle round bodied needle									
 <p>HR 65</p> <p>150 violet</p> <p>B0068919</p>									
 <p>HR 76</p> <p>100 violet</p> <p>B0068926</p>									
 <p>HR 76</p> <p>150 violet</p> <p>B0068927 B0068928</p>									

NOVOSYN®









ABSORBABLE SUTURES










Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 <p>HR 65</p>	240 violet					B0068929	B0068930		
 <p>HR 76</p>	150 violet						B0068931		
 <p>HR 65</p>	150 violet					B0068939	B0068940		
1/2 circle round bodied needle with short cutting point									
 <p>HRC 40</p>	90 violet				B0068950				
 <p>HRC 43s</p>	90 violet						B0068983		
 <p>HRC 48</p>	150 violet				B0068986	B0068987	B0068988		
 <p>HRC 40s</p>	90 violet						B0068992		

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 <p>HR 65</p>	70 undyed			B0069064					
1/2 circle reverse cutting needle									
 <p>HS 37</p>	90 undyed						B0069297		
 <p>HS 48</p>	90 undyed						B0069364		
1/2 circle round bodied needle with trocar point									
 <p>HRT 48</p>	150 undyed					B0069370			
1/2 circle round bodied needle									
 <p>HR 65</p>	140 undyed						B0069372		
 <p>HR 40s</p>	90 undyed						B0069558		
 <p>HR 37s</p>	90 undyed						B0069598		

NOVOSYN®










ABSORBABLE SUTURES









Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle with short cutting point										
 HRC 37s Easyslide	90 undyed				B0069886					
1/2 circle round bodied needle										
 HR 70 Easyslide	90 violet						B0088002			
 HR 70 Easyslide	4 x 70 90 violet						B0088003 B0088004			
1/2 circle round bodied needle with short cutting point										
 HRC 26 Easyslide	4 x 70 violet		B0088065							
1/2 circle round bodied needle										
 HR 48 Easyslide	5 x 45 violet			B0088350						
 HR 26 Easyslide	5 x 45 violet			B0088351						
 HR 37s Easyslide	5 x 45 violet				B0088352					
 HR 48 Easyslide	5 x 45 violet				B0088353	B0088354				

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 HR 40s take off	5 x 45 violet				B0088355				
 HR 22 take off	4 x 70 violet				B0088631				
 HR 37s take off	4 x 90 violet						B0088653		
 HR 65 take off	4 x 70 violet						B0088655	B0088656	
 HR 30 take off	4 x 70 violet				B0088661	B0088662			
 HR 40s take off	4 x 70 violet						B0088667	B0088668	
 HR 37 take off	4 x 70 violet				B0088670	B0088671			
 HR 65 take off	5 x 70 violet							B0088684	
 HR 48 take off	4 x 70 violet						B0088687	B0088688	

NOVOSYN®


ABSORBABLE SUTURES























Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 HR 26s take off	4 x 70 violet					B0088691	B0088692	B0088693	
 HR 37s take off	4 x 70 violet			B0088695			B0088697	B0088698	
1/2 circle reverse cutting needle									
 HS 37ss take off	3 x 45 violet					B0088701			
 HS 37s take off	4 x 45 violet						B0088705		
1/2 circle round bodied needle									
 HR 37ss take off	4 x 70 violet							B0088707	
 HR 40s take off	4 x 45 violet					B0088766	B0088767	B0088768	
 HR 43 take off	4 x 45 violet						B0088781	B0088782	
 HR 48 take off	4 x 45 violet			B0088785		B0088786	B0088787	B0088788	
 HR 26s take off	4 x 45 violet					B0088791	B0088792	B0088793	

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 HR 37s Easyslide take off	4 x 45 violet				B0088796	B0088797	B0088798		
 HR 30ss Easyslide take off	4 x 70 violet					B0088803			
 HR 37ss Easyslide take off	4 x 45 violet				B0088806	B0088807			
1/2 circle reverse cutting needle									
 HS 26ss Easyslide take off	4 x 45 violet				B0088906				
 HS 40ss Easyslide take off	3 x 45 4 x 45 violet					B0088911 B0088912			
1/2 circle round bodied needle									
 HR 37s Easyslide take off	4 x 70 undyed				B0089696				
 HR 48 Easyslide take off	4 x 45 undyed				B0089786				
 HR 37s Easyslide take off	4 x 45 undyed				B0089796	B0089797			

NOVOSYN®




















ABSORBABLE SUTURES












Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle reverse cutting needle									
	4 x 45 undyed					B0089912			
	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
Pre-cut sutures									
4 threads of 90 cm	4 x 90 violet				C0058064				
10 threads of 60 cm	10 x 60 violet	C0058180	C0058181	C0058182	C0058183				
1 thread of 140 cm	1 x 140 violet	C0058203	C0058204	C0058205	C0058207	C0058208	C0058209	C0058210	
2 threads of 70 cm	2 x 70 violet		C0058404	C0058405	C0058406	C0058407	C0058409	C0058410	
5 threads of 70 cm	5 x 70 violet		C0058444	C0058445	C0058446	C0058447			
5 threads of 90 cm	5 x 90 violet				C0058456	C0058457			
3 threads of 45 cm	3 x 45 violet	C0058623	C0058624	C0058625	C0058626	C0058627	C0058629	C0058630	
6 threads of 45 cm	6 x 45 violet	C0058643	C0058644	C0058645	C0058646				
12 threads of 45 cm	12 x 45 violet		C0058664	C0058665	C0058666				
1 thread of 140 cm	1 x 140 undyed			C0059205	C0059207	C0059208			
2 threads of 70 cm	2 x 70 undyed			C0059405			C0059409	C0059410	
5 threads of 70 cm	5 x 70 undyed			C0059445	C0059446	C0059447			
3 threads of 45 cm	3 x 45 undyed		C0059624	C0059625	C0059626	C0059627	C0059629	C0059630	
6 threads of 45 cm	6 x 45 undyed		C0059644	C0059645	C0059646				
12 threads of 45 cm	12 x 45 undyed		C0059664	C0059665	C0059666				

Needle length mm	Thread length cm colour	USP (metric)								
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	
1/2 circle round bodied needle										
 HR 10 	70 70 violet		C0068002	C0068003						
 HR 22 	20 violet							C0068004		
 HR 13 	70 violet		C0068006	C0068007	C0068008	C0068009				
 HR 17 	70 violet		C0068011	C0068012	C0068013	C0068014	C0068015	C0068016		
 HR 26 	20 violet							C0068018		
 HR 26s 	45 violet							C0068019		
 HR 6 	45 violet	C0068020								
 HR 37s 	45 violet								C0068021	
3/8 circle reverse cutting needle with precision point										
 DSMP 24 	45 violet							C0068022		
1/2 circle round bodied needle with trocar point										
 HRT 40 	70 violet							C0068023	C0068024	
3/8 circle reverse cutting needle with precision point										
 DSMP 24 	45 violet						C0068025			

NOVOSYN®















ABSORBABLE SUTURES



















Needle length mm	Thread length cm colour	USP (metric)								
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	
5/8 circle round bodied needle										
 FR 26b 	70 violet					C0068027				
	1/2 circle round bodied needle									
 HR 22 	70 violet			C0068028	C0068029					
		70 violet				C0068029N1				
		70 violet					C0068030			
		70 violet					C0068030N1			
		70 violet					C0068031			
		70 violet					C0068031N1			
		70 violet						C0068032		
 HR 22ss 	70 violet								C0068033	
Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
 HR 26 	70 violet							C0068034		
 HR 30 	70 violet		C0068035							
 HR 22s 	70 violet				C0068036					
		70 violet				C0068036N1				
		70 violet					C0068037			
		70 violet					C0068037N1			
		70 violet					C0068038			
		70 violet					C0068038N1			

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
 HR 26	70 violet	C0068039	C0068040						
	 70 violet		C0068040N1						
	70 violet			C0068041					
	 70 violet			C0068041N1					
	70 violet				C0068042				
	 70 violet				C0068042N1				
	70 violet					C0068043			
	 70 violet					C0068043N1			
	70 violet						C0068044		
	 70 violet						C0068044N1		
 HR 30	70 violet			C0068046				C0068045	
	 70 violet			C0068046N1					
	70 violet				C0068047				
	 70 violet				C0068047N1				
	70 violet					C0068048			
	 70 violet					C0068048N1			
	70 violet						C0068049		
	 70 violet						C0068049N1		

NOVOSYN®










ABSORBABLE SUTURES













Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle										
 <p>HR 37</p>	70 violet			C0068050						
	 70 violet			C0068050N1						
	70 violet				C0068051					
	 70 violet				C0068051N1					
	70 violet					C0068052				
	 70 violet					C0068052N1				
	70 violet						C0068053			
	 70 violet						C0068053N1			
	 70 violet							C0068054		
	 <p>HR 40s</p>	70 violet				C0068055				
 70 violet					C0068055N1					
70 violet						C0068056				
 70 violet						C0068056N1				
70 violet							C0068057			
 70 violet							C0068057N1			
 70 violet				C0068059				C0068058		
 <p>HR 48</p>		70 violet				C0068060				
		 70 violet				C0068060N1				
		 70 violet					C0068061			

Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle										
 HR 48 	70 violet					C0068061N1				
	70 violet							C0068062		
	70 violet							C0068062N1		
	70 violet								C0068063	
 HR 60 	70 violet				C0068064	C0068065	C0068066			
 HR 26 	120 violet				C0068067					
 HR 60 	70 violet								C0068068	
 HR 22 	20 violet			C0068069						
 HR 43s 	70 violet								C0068070	
 HR 22 	120 violet			C0068078						
 HR 26 	20 violet			C0068079						
1/2 circle round bodied needle with short cutting point										
 HRC 17 	70 violet							C0068081		

NOVOSYN®







ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle										
 HR 48	70 violet									
	70 violet					C0068061N1				
	70 violet						C0068062			
	70 violet						C0068062N1			
	70 violet							C0068063		
	70 violet									
 HR 60	70 violet					C0068064	C0068065	C0068066		
 HR 26	120 violet					C0068067				
 HR 60	70 violet								C0068068	
 HR 22	20 violet				C0068069					
 HR 43s	70 violet								C0068070	
 HR 22	120 violet				C0068078					
 HR 26	20 violet				C0068079					
1/2 circle round bodied needle with short cutting point										
 HRC 17	70 violet							C0068081		

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle reverse cutting needle									
 HS 26ss Easyslide	70 violet								C0068084
 HS 37ss Easyslide	70 violet				C0068085	C0068086	C0068087		
1/2 circle round bodied needle									
 HR 26 Easyslide	120 violet			C0068088					
 HR 26s Easyslide	70 violet			C0068089					
	 70 violet			C0068089N1					
	 70 violet				C0068090				
	 70 violet				C0068090N1				
	 70 violet					C0068091			
	 70 violet					C0068091N1			
	 70 violet						C0068092		
	 70 violet						C0068092N1		
	 70 violet							C0068093	

Needle length mm	Thread length cm colour	USP (metric)								
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	







1/2 circle round bodied needle

 	45 violet							C0068109		
 	70 violet							C0068110	C0068111	
 	70 violet				C0068113	C0068114	C0068115	C0068116		

3/8 circle round bodied needle









 	70 violet						C0068117			
---	--------------	--	--	--	--	--	----------	--	--	--

1/2 circle round bodied needle

 	70 violet									C0068118
 	70 violet					C0068120				
 	70 violet							C0068121	C0068122	C0068123









Needle length mm	Thread length cm colour	USP (metric)								
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	

3/8 circle round bodied needle

 	30 violet	C0068133								
 	45 violet				C0068134					
 	45 violet					C0068137				
 	45 violet					C0068138				


NOVOSYN[®]

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
1/2 circle round bodied needle with short cutting point									
 HRC 37	90 violet						C0068139	C0068140	
	90 violet							C0068140N1	
	90 violet								C0068141
	90 violet								C0068141N1
1/2 circle round bodied needle									
 HRC 37	90 violet								C0068142
	90 violet								C0068142N1
1/2 circle round bodied needle									
 HR 17b	70 violet			C0068144	C0068147				
	 HR 6	30 violet	C0068149						
3/8 circle round bodied needle									
 DR 18	70 violet				C0068150	C0068151			
1/2 circle round bodied needle									
 HR 17b	70 violet					C0068152			
Ski Needle round bodied needle									
 SKR 19	45 violet					C0068154			
3/8 circle round bodied needle									
 DR 18	70 violet			C0068155					

Needle length mm	Thread length cm colour	USP (metric)								
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	


3/8 circle round bodied needle

 DR 22	70 violet				C0068159	C0068160			
	70 violet					C0068160N1			
	70 violet						C0068161		


3/8 circle reverse cutting needle with precision point

 2xDSMP 13	45 violet			C0068162					
--	--------------	--	--	----------	--	--	--	--	--

1/2 circle round bodied needle

 2xHR 17b	70 violet			C0068163					
---	--------------	--	--	----------	--	--	--	--	--

3/8 circle round bodied needle

 DR 26	70 violet					C0068164	C0068165		
	70 violet						C0068165N1		


Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	

 DR 32	70 violet								C0068168
--	--------------	--	--	--	--	--	--	--	----------

1/2 circle round bodied needle










 HR 43s	70 violet							C0068169	
---	--------------	--	--	--	--	--	--	----------	--





















3/8 circle round bodied needle















 DR 32	70 violet				C0068170	C0068171			
	70 violet					C0068171N1			
	70 violet						C0068172	C0068173	















NOVOSYN[®]








ABSORBABLE SUTURES

















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 40	70 violet						C0068175			
	70 violet						C0068175N1			
 2xHR 17b	70 violet			C0068176						
3/8 circle round bodied needle										
 DR 26	120 violet				C0068179					
1/2 circle round bodied needle										
 HR 22b	70 violet			C0068181						
Straight round bodied needle										
 GR 51	90 violet					C0068185				
1/2 circle round bodied needle										
 HR 22b	70 violet				C0068186					
	70 violet				C0068186N1					
1/2 circle round bodied needle with short cutting point										
 HRC 37	150 violet					C0068187				
1/2 circle round bodied needle										
 HR 22b	70 violet					C0068189				
 2xHR 22b	70			C0068190						
	90			C0068191						
	70 violet				C0068192					

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Hook round bodied needle with short cutting point										
 JRC 30s 	120 violet						C0068193	C0068194		
1/2 circle round bodied needle										
 HR 26b  	70 violet			C0068196	C0068198					
	70 violet				C0068198N1					
3/8 circle reverse cutting needle										
 DS 9 	45 violet	C0068202								
1/2 circle reverse cutting needle										
 HS 15 	45 violet		C0068203	C0068204						
 HS 18 	70 violet							C0068205		
3/8 circle reverse cutting needle										
 DS 12 	45 violet	C0068206	C0068207	C0068208						
1/2 circle round bodied needle										
 HR 26b  	70 violet					C0068209				
	70 violet					C0068209N1				
 2xHR 26b 	70 violet				C0068210					
Straight reverse cutting needle										
 GS 13 	45 violet			C0068211						

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
3/8 circle reverse cutting needle										
 DS 24 	70 violet		C0068232							
1/2 circle reverse cutting needle										
 HS 15 	70 violet		C0068233							
3/8 circle reverse cutting needle										
 DS 24 	70 violet			C0068234						
		70 violet		C0068234N1						
		70 violet			C0068235					
		70 violet			C0068235N1					
		70 violet				C0068236				
		70 violet				C0068236N1				
		70 violet					C0068237			
		70 violet						C0068237N1		
		90 violet						C0068239	C0068240	
1/2 circle reverse cutting needle										
 HS 26 	70 violet				C0068241					
		70 violet				C0068241N1				
		70 violet					C0068242			
		70 violet						C0068242N1		
		70 violet							C0068243	C0068244









Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 HS 40s 	70 violet						C0068267	C0068268		
	70 violet							C0068268N1		
	70 violet								C0068269	
	70 violet								C0068269N1	
Straight reverse cutting needle										
 GS 51 	70 violet				C0068272	C0068274	C0068275			
Straight round bodied needle										
 GR 19 	20 violet			C0068281						
3/8 circle reverse cutting needle										
 DS 39 	70 violet						C0068282	C0068283		
Straight round bodied needle										
 GR 19 	20 70 violet				C0068284 C0068285					
Ski Needle round bodied needle										
 SKR 26b 	20 violet					C0068287				
Straight round bodied needle										
 GR 19b 	20 violet			C0068288						

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
3/8 circle round bodied needle with short cutting point										
 DRC 26 Easyslide	70 violet						C0068317			
3/8 circle round bodied needle with trocar point										
 DRT 18 Easyslide	70 violet		C0068319	C0068320	C0068321					
3/8 circle round bodied needle with short cutting point										
 DRC 32 Easyslide	70 70 violet						C0068322	C0068324		
1/2 circle round bodied needle with short cutting point										
 HRC 22 Easyslide	70 violet			C0068325	C0068326					
 HRC 37ss Easyslide	70 violet						C0068329			
 HRC 22 Easyslide	70 violet						C0068331	C0068332	C0068333	
 HRC 26 Easyslide	20 70 violet				C0068334 C0068335					
	70 violet				C0068335N1					
	70 violet					C0068336				
	70 violet					C0068336N1				
	20 violet					C0068337				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
 HRC 40ss 	70 violet									C0068359
	70 violet				C0068360	C0068361				
 HRC 37 	70 violet									
	70 violet					C0068361N1				
	70 violet						C0068362	C0068363		
1/2 circle round bodied needle										
 HR 22b 	120 violet				C0068365					
1/2 circle round bodied needle with short cutting point										
 HRC 43ss 	70 violet									C0068368
	70 violet						C0068372	C0068371		
 HRC 30 	70 violet									
	70 violet						C0068372N1			
	70 violet				C0068373				C0068374	
 HRC 30s 	90 violet									C0068375
3/8 circle round bodied needle with trocar point										
 DRT 22 	70 violet			C0068376	C0068377					
3/8 circle reverse cutting needle										
 DS 45 	70 violet						C0068379			

NOVOSYN®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 40 Easyslide	70 violet						C0068381			
3/8 circle round bodied needle with trocar point										
 DRT 26 Easyslide	70 violet				C0068382	C0068383				
1/2 circle round bodied needle with trocar point										
 HRT 37s Easyslide	90 violet				C0068384					
3/8 circle round bodied needle with trocar point										
 DRT 32 Easyslide	70 violet					C0068386	C0068387			
1/2 circle round bodied needle with trocar point										
 HRT 37 Easyslide	70 violet					C0068392	C0068393	C0068391		
 HRT 37s Easyslide	70 violet					C0068395	C0068396	C0068397		
5/8 circle reverse cutting needle										
 FS 25 Easyslide	70 violet				C0068398					
1/2 circle round bodied needle with trocar point										
 HRT 48 Easyslide	70 violet					C0068400	C0068401	C0068402	C0068403	

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)

Straight round bodied needle



GR 51




70
violet

C0068405

1/2 circle round bodied needle with short cutting point



HRC 43ss



70
violet


C0068406

Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)

3/8 circle reverse cutting needle with precision point



DSMP 7



70
violet

C0068407

Hook round bodied needle



JR 36ss



70
violet

C0068408

3/8 circle round bodied needle



DR 40



70
violet

C0068409

Hook round bodied needle



JR 36ss



70
violet

C0068410

1/2 circle round bodied needle with trocar point



HRT 40s



































70
violet

















C0068413

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

NOVOSYN®


















ABSORBABLE SUTURES









Needle length mm	Thread length cm colour	USP (metric)																																																																																																														
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)																																																																																																							
5/8 circle round bodied needle																																																																																																																
 FR 26 	70 violet					C0068414	C0068415	C0068416	C0068417																																																																																																							
	3/8 circle reverse cutting needle																																																																																																															
 DS 19 	70 violet			C0068419	C0068420																																																																																																											
	 70 violet				C0068420N1																																																																																																											
	70 violet					C0068421																																																																																																										
	 70 violet					C0068421N1																																																																																																										
	70 violet						C0068422																																																																																																									
	 70 violet						C0068422N1																																																																																																									
5/8 circle round bodied needle																																																																																																																
 FR 36 	70 violet					C0068424	C0068425	C0068426	C0068427																																																																																																							
	 FR 26s 	70 violet								C0068428																																																																																																						
<table border="1"> <thead> <tr> <th rowspan="2">Needle length mm</th> <th rowspan="2">Thread length cm colour</th> <th colspan="8">USP (metric)</th> </tr> <tr> <th>6/0 (0.7)</th> <th>5/0 (1)</th> <th>4/0 (1.5)</th> <th>3/0 (2)</th> <th>2/0 (3)</th> <th>0 (3.5)</th> <th>1 (4)</th> <th>2 (5)</th> </tr> </thead> <tbody> <tr> <td colspan="10">1/2 circle reverse cutting needle</td> </tr> <tr> <td rowspan="6">  HS 23  </td> <td>70 violet</td> <td></td> <td></td> <td></td> <td>C0068430</td> <td></td> <td></td> <td></td> <td>C0068429</td> </tr> <tr> <td> 70 violet</td> <td></td> <td></td> <td></td> <td>C0068430N1</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>70 violet</td> <td></td> <td></td> <td></td> <td></td> <td>C0068431</td> <td></td> <td></td> <td></td> </tr> <tr> <td> 70 violet</td> <td></td> <td></td> <td></td> <td></td> <td>C0068431N1</td> <td></td> <td></td> <td></td> </tr> <tr> <td>70 violet</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>C0068432</td> <td></td> <td></td> </tr> <tr> <td> 70 violet</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td colspan="10">3/8 circle round bodied needle</td> </tr> <tr> <td>  DR 8  </td> <td>30 violet</td> <td>C0068433</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>										Needle length mm	Thread length cm colour	USP (metric)								6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	1/2 circle reverse cutting needle										 HS 23 	70 violet				C0068430				C0068429	 70 violet				C0068430N1					70 violet					C0068431				 70 violet					C0068431N1				70 violet						C0068432			 70 violet									3/8 circle round bodied needle										 DR 8 	30 violet	C0068433							
Needle length mm	Thread length cm colour	USP (metric)																																																																																																														
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)																																																																																																							
1/2 circle reverse cutting needle																																																																																																																
 HS 23 	70 violet				C0068430				C0068429																																																																																																							
	 70 violet				C0068430N1																																																																																																											
	70 violet					C0068431																																																																																																										
	 70 violet					C0068431N1																																																																																																										
	70 violet						C0068432																																																																																																									
	 70 violet																																																																																																															
3/8 circle round bodied needle																																																																																																																
 DR 8 	30 violet	C0068433																																																																																																														

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
3/8 circle round bodied needle										
 DR 12 	70 violet	C0068434	C0068435	C0068436						
1/2 circle round bodied needle										
 HR 22 	90 violet			C0068438						
1/2 circle reverse cutting needle										
 HS 18 	70 violet		C0068439							
 HS 26s 	90 violet					C0068440				
1/2 circle round bodied needle with short cutting point										
 HRC 26s 	70 violet						C0068441			
1/2 circle reverse cutting needle										
 HS 26 	90 violet					C0068442				
1/2 circle round bodied needle										
 HR 22 	70 violet								C0068443	
1/2 circle round bodied needle with short cutting point										
 HRC 26s 	70 violet							C0068444		

NOVOSYN®








ABSORBABLE SUTURES













Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 30  	90 violet				C0068446	C0068447			
	90 violet					C0068447N1			
	90 violet						C0068448		
	90 violet						C0068448N1		
 HR 13 	45 90 violet	C0068449	C0068450	C0068451 C0068452					
3/8 circle reverse cutting needle									
 DS 48 	70 violet							C0068453	
1/2 circle round bodied needle									
 HR 40  	70 violet							C0068457	
	70 violet							C0068457N1	
1/2 circle reverse cutting needle									
 HS 26 	90 violet				C0068458				
3/8 circle round bodied needle with short cutting point									
 DRC 18 	70 violet			C0068460					
1/2 circle reverse cutting needle									
 HS 48  	90 violet						C0068461	C0068462	
	90 violet							C0068462N1	

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 HS 48f	90 violet									C0068464
Straight reverse cutting needle										
 GS 60	45 violet				C0068465					
1/2 circle round bodied needle										
 HR 40s	90 violet				C0068466	C0068467				
	90 violet					C0068467N1				
3/8 circle round bodied needle with short cutting point										
 DRC 18	70 violet				C0068468					
1/2 circle round bodied needle with trocar point										
 HRT 43s	90 violet									C0068470
 HRT 43	90 violet						C0068471	C0068472		
 HRT 43s	70 violet									C0068474
Hook round bodied needle										
 JR 30	70 violet						C0068475			

NOVOSYN[®]



















ABSORBABLE SUTURES


















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Hook round bodied needle with short cutting point										
 JRC 30s	70 violet						C0068477	C0068476		
		Hook round bodied needle								
 JR 40s	70 violet							C0068478		
		Hook round bodied needle with short cutting point								
 2xJRC 30s	45 violet						C0068479			
		1/2 circle reverse cutting needle								
 HS 26s	90 violet						C0068480			
		Hook round bodied needle								
 JR 30s	70 violet							C0068481		
		Hook round bodied needle with short cutting point								
 2xJRC 30s	60 violet					C0068483				
		1/2 circle round bodied needle with trocar point								
 HRT 40s	70 violet							C0068484		

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 HS 40 	90 violet						C0068485	C0068486		
	90 violet							C0068486N1		
	90 violet								C0068487	
	90 violet								C0068487N1	
5/8 circle round bodied needle										
 FR 40 	70 violet						C0068489	C0068490	C0068491	
1/2 circle round bodied needle with trocar point										
 HRT 40s 	70 violet									C0068492
1/2 circle reverse cutting needle										
 HS 30 	90 violet								C0068493	
 HS 37 	90 violet				C0068494	C0068495				
	90 violet						C0068495N1			
	90 violet							C0068496		
	90 violet							C0068496N1		
	90 violet								C0068497	
90 violet									C0068497N1	
1/2 circle round bodied needle with trocar point										
 HRT 48 	90 violet				C0068504	C0068500	C0068501	C0068502		

NOVOSYN®









ABSORBABLE SUTURES







Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
5/8 circle round bodied needle with short cutting point										
 	90 violet					C0068505				
3/8 circle reverse cutting needle										
 	70 violet	C0068506	C0068507	C0068508						
 	90 70 violet		C0068511 C0068512	C0068513	C0068514					
1/2 circle round bodied needle with protection point										
 	90 violet								C0068516	
1/2 circle round bodied needle with blunt point										
 	70 violet					C0068517	C0068518			
1/2 circle reverse cutting needle										
 	70 violet			C0068522						
5/8 circle round bodied needle with short cutting point										
 	90 violet			C0068523						
1/2 circle reverse cutting needle										
 	90 violet					C0068524				
1/2 circle round bodied needle with trocar point										
 	70 violet								C0068525	

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 22b 	90 violet					C0068526				
3/8 circle reverse cutting needle with precision point										
 DSMP 19 	70 violet		C0068528							
 DSMP 16 	45 violet	C0068530		C0068531	C0068529					
1/2 circle round bodied needle with trocar point										
 HRT 37s 	90 violet						C0068532	C0068533	C0068534	
3/8 circle reverse cutting needle										
 DS 24 	45 violet			C0068536						
3/8 circle reverse cutting needle with precision point										
 DSMP 16 	45 violet		C0068538							
1/2 circle round bodied needle										
	90 violet			C0068540	C0068541					
	 90 violet				C0068541N1					
	90 violet					C0068542				
	 90 violet					C0068542N1				
	90 violet						C0068543			
	 90 violet						C0068543N1			
	90 violet							C0068544	C0068545	
	 90 violet									

NOVOSYN®






ABSORBABLE SUTURES


















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 30s 	70 violet						C0068547	C0068548		
	70 violet							C0068548N1		
	70 violet								C0068549	
	70 violet								C0068549N1	
	90 violet					C0068550				
	90 violet					C0068550N1				
 HR 37 	90 violet					C0068551				
	90 violet					C0068551N1				
	90 violet						C0068552			
	90 violet						C0068552N1			
	90 violet							C0068553		
	90 violet							C0068553N1		
 HR 40s 	90 violet						C0068556			
	90 violet						C0068556N1			
	90 violet							C0068557		
	90 violet							C0068557N1		
 HR 30ss 	70 violet							C0068559		

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 48	90 violet						C0068560			
	90 violet						C0068560N1			
	90 violet							C0068561		
	90 violet							C0068561N1		
	90 violet								C0068562	
	90 violet									C0068562N1
3/8 circle round bodied needle										
 DR 32	90 violet						C0068566	C0068567	C0068568	
1/2 circle round bodied needle										
 HR 40	90 violet							C0068569		
	90 violet							C0068569N1		
1/2 circle round bodied needle with short cutting point										
 HRC 48	90 violet					C0068570				
1/2 circle round bodied needle										
 HR 40	90 violet								C0068571	
	90 violet									C0068571N1
1/2 circle round bodied needle with short cutting point										
 HRC 48	90 violet						C0068572	C0068573	C0068574	

NOVOSYN[®]

















ABSORBABLE SUTURES









Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 37ss	90 violet						C0068576	C0068577		
	90 violet								C0068577N1	
 HR 26ss	90 violet							C0068580		
1/2 circle reverse cutting needle										
 HS 37s	90 violet						C0068582			
	90 violet						C0068582N1			
	90 violet							C0068583		
	90 violet							C0068583N1		
1/2 circle round bodied needle										
 HR 37s	150 violet						C0068584			
1/2 circle reverse cutting needle										
 HS 40s	90 violet						C0068585	C0068586		
	90 violet							C0068586N1		
	90 violet								C0068587	
	90 violet								C0068587N1	

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 26s	90 violet				C0068589	C0068590			
	 90 violet					C0068590N1			
	90 violet						C0068591		
	 90 violet						C0068591N1		
	90 violet							C0068592	
	 90 violet							C0068592N1	
1/2 circle reverse cutting needle									
 HS 37s	90 violet				C0068594				
									
1/2 circle round bodied needle									
 HR 37s	90 violet					C0068595			
	  90 violet						C0068595N1		
1/2 circle reverse cutting needle									
 HS 37s	90 violet					C0068596			
	  90 violet						C0068596N1		
1/2 circle round bodied needle									
 HR 37s	90 violet							C0068597	
	  90 violet							C0068597N1	
1/2 circle reverse cutting needle									
 HS 26	120 violet				C0068600				
									

NOVOSYN[®]









ABSORBABLE SUTURES

















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 HS 26s 	120									
	violet					C0068603				
1/2 circle round bodied needle with blunt point										
 HRN 30 	70									
	violet					C0068604				
 HRN 35 	70									
	violet								C0068608	
1/2 circle reverse cutting needle										
 HS 26s 	120									
	violet								C0068610	
1/2 circle round bodied needle										
 HR 26s 	90									
	violet									C0068613
1/2 circle round bodied needle with blunt point										
 HRN 35 	90									
	violet						C0068617	C0068618	C0068619	
3/8 circle reverse cutting needle with precision point										
 DSMP 13 	45									
	violet	C0068622	C0068623	C0068624						
1/2 circle round bodied needle with short cutting point										
 HRC 43s 	70									
	violet							C0068625		

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
3/8 circle reverse cutting needle with precision point										
 DSMP 19 Easyslide	45 violet		C0068626							
1/2 circle round bodied needle with short cutting point										
 HRC 43s Easyslide	70 violet								C0068627	C0068628
Ski Needle round bodied needle										
 SKR 26 Easyslide	30 violet						C0068629			
3/8 circle reverse cutting needle with precision point										
 DSMP 19 Easyslide	45 violet			C0068630	C0068631					
1/2 circle round bodied needle										
 HR 30s Easyslide	90 violet								C0068633	
1/2 circle round bodied needle with blunt point										
 HRN 40 Easyslide	70 90 violet							C0068634 C0068636	C0068635 C0068637	
 HRN 30 Easyslide	70 violet				C0068638					
1/2 circle reverse cutting needle										
 HS 26s Easyslide	70 violet				C0068641					

NOVOSYN®












ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 26ss Easyslide	90 violet									C0068642
 HR 30 Easyslide	90 violet								C0068645	
 HR 30ss Easyslide	90 violet								C0068646	
1/2 circle round bodied needle with blunt point										
 HRN 45 Easyslide	90 violet								C0068647	
1/2 circle round bodied needle with trocar point										
 HRT 37 Easyslide	90 violet				C0068650	C0068651	C0068652	C0068653		
3/8 circle reverse cutting needle with precision point										
 DSMP 39 Easyslide	70 violet				C0068659	C0068660				
1/2 circle round bodied needle with blunt point										
 HRN 50 Easyslide	90 violet							C0068663	C0068664	
3/8 circle reverse cutting needle with precision point										
 DSMP 30 Easyslide	70 violet				C0068669	C0068670	C0068671			

Needle length mm	Thread length cm colour	USP (metric)								
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	0 (3.5)
1/2 circle reverse cutting needle with precision point										
 2xHSMP 9 	35 violet	C0068673								
3/8 circle reverse cutting needle with precision point										
 DSMP 11 	45 70 violet		C0068684	C0068685 C0068688						
 DSMP 35 	70 violet							C0068690		
1/2 circle reverse cutting needle										
 HS 48 	70 violet							C0068691		
1/2 circle round bodied needle with short cutting point										
 HRC 26 	70 violet					C0068694				
Needle length mm	Thread length cm colour	USP (metric)								
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	
Straight round bodied needle										
 GR 19 	70 violet									C0068695
3/8 circle reverse cutting needle										
 2xDS 12 	45 violet		C0068706							
3/8 circle reverse cutting needle with precision point										
 2xDSMP 7 	45 violet		C0068708							



NOVOSYN®

ABSORBABLE SUTURES



Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle with short cutting point									
 Easyslide	70 violet								C0068716
 Easyslide	70 violet								C0068721
3/8 circle reverse cutting needle with precision point									
 Easyslide	70 violet			C0068723					
 Easyslide	45 violet		C0068726						
1/4 circle lancet needle with micro-point									
 Easyslide	45 violet	C0068732							
Ski Needle round bodied needle									
 Easyslide	120 violet					C0068744			
 Easyslide	40 violet					C0068745			
 Easyslide	20 violet					C0068746			
1/2 circle round bodied needle									
 Easyslide	70 violet							C0068747	
Ski Needle round bodied needle									
 Easyslide	20 120 violet							C0068749 C0068751	
 Easyslide	20 120 violet				C0068752 C0068755	C0068753 C0068756	C0068754 C0068757		

Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)


1/2 circle reverse cutting needle with precision point

 2xHSMP 9 	45								
	violet		C0068758						

1/2 circle round bodied needle

  HR 20	90					C0068763			
	70				C0068765				
	violet								
	70				C0068765N1				
	violet								
	70					C0068766			
violet									
70						C0068766N1			
violet									
70							C0068767		
violet									

3/8 circle reverse cutting needle

 2xDS 12 	45								
	violet				C0068769				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)

1/2 circle reverse cutting needle

  HS 40s	70								
	violet								C0068770

5/8 circle round bodied needle

  FR 26	90					C0068775			
	violet								

Ski Needle round bodied needle





















  SKR 26	120						C0068781		
	violet								

















3/8 circle reverse cutting needle with precision point

  DSMP 16	70								
	violet		C0068782	C0068783	C0068784				

NOVOSYN®



















ABSORBABLE SUTURES

















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
 	70 violet						C0068785			
1/2 circle round bodied needle										
 	70 violet							C0068791		C0068793
5/8 circle round bodied needle with short cutting point										
 	30 violet	C0068798								
1/2 circle reverse cutting needle										
 	70 violet							C0068805	C0068806	
Progressive curved round bodied needle with short cutting point										
 	120 violet						C0068807			
5/8 circle round bodied needle										
 	90 violet							C0068810		
 	70 90 violet			C0068813	C0068818	C0068815 C0068820	C0068817			
3/8 circle reverse cutting needle with precision point										
 	70 violet			C0068823	C0068824					
 	70 violet					C0068824N1				
3/8 circle round bodied needle with short cutting point										
 	70 violet					C0068827				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
5/8 circle round bodied needle										
 	90 violet							C0068828		
3/8 circle round bodied needle with short cutting point										
 	20 violet						C0068830			
 	20 violet				C0068833					
1/2 circle reverse cutting needle										
 	70 violet									C0068840
1/2 circle round bodied needle with short cutting point										
 	70 violet						C0068842	C0068843	C0068844	
1/2 circle round bodied needle										
 	90 violet									C0068847
	90 violet									C0068847N1
1/2 circle round bodied needle with short cutting point										
 	140 violet					C0068852				
1/2 circle round bodied needle										
 	90 violet			C0068853	C0068854	C0068855				

NOVOSYN[®]



















ABSORBABLE SUTURES








Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Hook round bodied needle with short cutting point										
 JRC 30s 	90 violet									C0068857
5/8 circle round bodied needle										
 FR 36 	90 violet								C0068860	
 2xFR 36 	90 violet								C0068861	
3/8 circle reverse cutting needle with precision point										
 DSMP 24  	70 violet				C0068864	C0068865				
	70 violet						C0068865N1			
5/8 circle round bodied needle										
 FR 40 	90 violet							C0068866	C0068867	
1/2 circle round bodied needle										
 2xHR 17 	90 violet		C0068869	C0068870						
 2xHR 22 	90 violet				C0068874	C0068875				
 HR 22  	90 violet				C0068877					
	90 violet						C0068877N1			
	120 90 violet							C0068880 C0068881		

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 17f 	70 violet			C0068888						
1/2 circle round bodied needle with short cutting point										
 HRC 37s 	90 violet				C0068893					
	 90 violet				C0068893N1					
	90 violet					C0068894				
	 90 violet					C0068894N1				
	90 violet						C0068895			
	 90 violet						C0068895N1			
	90 violet							C0068896		
	 90 violet							C0068896N1		
1/2 circle round bodied needle										
 2xHR 26 	90 violet				C0068901	C0068902				
 2xHR 17 	70 violet				C0068904					
 2xHR 26 	70 violet						C0068907			
1/2 circle round bodied needle with short cutting point										
 HRC 40 	70 violet					C0068908	C0068909			

NOVOSYN®








ABSORBABLE SUTURES









Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Straight reverse cutting needle										
 2xGS 13 	45 violet			C0068920						
1/2 circle round bodied needle										
 HR 22s 	70 violet			C0068932	C0068934					
 2xHR 22s 	70 violet						C0068937			
Straight round bodied needle										
 2xGR 19 	35 violet			C0068943						
 GR 19 	35 violet			C0068944						
1/2 circle round bodied needle										
 HR 43 	70 violet						C0068946	C0068947		
	70 violet							C0068947N1		
	70 violet									C0068948
1/2 circle round bodied needle with short cutting point										
 HRC 40 	90 violet					C0068949		C0068951		
1/2 circle round bodied needle										
 HR 37s 	90 violet			C0068953	C0068955		C0068956			
	90 violet						C0068956N1			
1/2 circle reverse cutting needle										
 HS 26 	120 violet					C0068957				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
 HRC 40s Easyslide	70 violet						C0068958	C0068959	C0068960	
	70 violet								C0068960N1	
	70 violet									C0068961
 HRC 48 Easyslide	70 violet						C0068962	C0068963	C0068964	C0068965
 HRC 43 Easyslide	70 violet						C0068971	C0068972		
1/2 circle reverse cutting needle										
 HS 26ss Easyslide	70 violet							C0068975		
1/2 circle round bodied needle with short cutting point										
 HRC 43 Easyslide	90 violet						C0068976	C0068977	C0068978	
 HRC 43s Easyslide	90 violet								C0068982	
Straight round bodied needle										
 GR 65 Easyslide	90 violet						C0068994			

NOVOSYN®










ABSORBABLE SUTURES




















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with trocar point										
 HRT 43 ⊕ Easyslide	90 violet					C0068997				
5/8 circle round bodied needle with short cutting point										
 FR 17 ⊙ Easyslide	70 undyed			C0069003						
1/2 circle round bodied needle with short cutting point										
 HRC 37s ⊖ Easyslide	90 undyed					C0069004				
1/2 circle round bodied needle										
 HR 13 ⊙ Easyslide	70 undyed		C0069007	C0069008						
1/2 circle round bodied needle with short cutting point										
 HRC 40ss ⊖ Easyslide	70 undyed					C0069010				
1/2 circle round bodied needle										
 HR 17 ⊙ Easyslide	70 undyed	C0069011	C0069012	C0069013	C0069014	C0069015				
 HR 22 ⊙ Easyslide	70 undyed			C0069029	C0069030	C0069031		C0069033		

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 26	70 undyed			C0069040	C0069041					
	70 undyed				C0069041N1					
	70 undyed					C0069042				
	70 undyed					C0069042N1				
	70 undyed						C0069043			
 HR 30	70 undyed					C0069047	C0069048			
 HR 37	70 undyed				C0069050	C0069051	C0069052			
 HR 40s	70 undyed					C0069055		C0069057		
 HR 48	70 undyed				C0069059					
1/2 circle reverse cutting needle										
 HS 26ss	70 undyed							C0069083		
 HS 26s	70 undyed						C0069084			
 HS 37ss	70 undyed					C0069085	C0069088	C0069087		

NOVOSYN[®]























ABSORBABLE SUTURES



















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 26s <small>Easyslide</small>	70 undyed					C0069089	C0069090	C0069091	C0069092	
1/2 circle reverse cutting needle										
 HS 40s <small>Easyslide</small>	70 undyed								C0069093	
1/2 circle round bodied needle										
 HR 37s <small>Easyslide</small>	70 undyed					C0069099	C0069095	C0069096	C0069097	C0069098
Progressive curved round bodied needle										
 PR 16 <small>Easyslide</small>	45 undyed			C0069135						
3/8 circle reverse cutting needle										
 DS 30 <small>Easyslide</small>	90 undyed							C0069148		
1/2 circle reverse cutting needle										
 HS 37s <small>Easyslide</small>	90 70 undyed						C0069187 C0069195	C0069196	C0069197	
 HS 15 <small>Easyslide</small>	45 undyed		C0069203	C0069204						
3/8 circle reverse cutting needle										
 DS 12 <small>Easyslide</small>	45 undyed	C0069206	C0069207	C0069208						
 DS 16 <small>Easyslide</small>	45 undyed		C0069212	C0069213	C0069214					

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
3/8 circle reverse cutting needle										
 DS 19f 	45 undyed			C0069217						
 DS 19 	45 undyed		C0069219	C0069220						
	45 undyed			C0069220N1						
	45 undyed				C0069221					
	45 undyed				C0069221N1					
	90 undyed			C0069223						
1/2 circle reverse cutting needle										
 HS 23 	45 undyed				C0069230					
3/8 circle reverse cutting needle										
 DS 24 	90 70 undyed			C0069234	C0069235	C0069232				
	70 undyed				C0069235N1					
	70 undyed					C0069236				
	70 undyed					C0069236N1				
	70 undyed						C0069237			
 DS 30 	70 undyed				C0069246	C0069247		C0069244		
	70 undyed					C0069247N1				
	70 undyed						C0069248			
 DS 35 	70 undyed					C0069256				
Straight reverse cutting needle										
 GS 60 	70 undyed				C0069265					

NOVOSYN®









ABSORBABLE SUTURES



















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Straight reverse cutting needle										
 GS 60f 	70 undyed					C0069266				
 GS 60 	70 undyed						C0069267			
 GS 51 	70 undyed			C0069271	C0069272	C0069273				
 GS 60 	70 undyed			C0069280						
3/8 circle reverse cutting needle										
 DS 48 	70 undyed									C0069282
1/2 circle reverse cutting needle										
 HS 26s 	70 undyed						C0069290		C0069293	
 HS 37 	90 undyed						C0069294	C0069295	C0069296	
5/8 circle round bodied needle with short cutting point										
 FRC 22 	90 undyed			C0069299	C0069298					
3/8 circle reverse cutting needle										
 DS 16 	70 undyed					C0069302				
1/2 circle round bodied needle with short cutting point										
 HRC 17 	45 undyed		C0069309							
3/8 circle round bodied needle with trocar point										
 DRT 18 	70 undyed			C0069320	C0069321					

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 HS 18 	70 undyed				C0069322					
1/2 circle round bodied needle with short cutting point										
 HRC 22 	70 undyed				C0069324					
3/8 circle reverse cutting needle with precision point										
 DSMP 16 	70 undyed		C0069325	C0069326	C0069327					
3/8 circle reverse cutting needle										
 DS 12 	70 undyed	C0069330	C0069331							
1/2 circle round bodied needle with short cutting point										
 HRC 22 	70 undyed				C0069335					
 HRC 17 	70 undyed		C0069339	C0069340	C0069341					
1/2 circle reverse cutting needle										
 HS 48 	90 undyed						C0069362	C0069363		
1/2 circle round bodied needle with trocar point										
 HRT 43s 	90 undyed						C0069378			
 HRT 37 	70 undyed						C0069392			

NOVOSYN®























ABSORBABLE SUTURES










Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with trocar point										
 HRT 40 ⊕ Easyslide	70 undyed					C0069400				
1/2 circle round bodied needle										
 HR 13 ● Easyslide	45 undyed	C0069405	C0069407	C0069408						
5/8 circle round bodied needle										
 FR 26 ● Easyslide	70 undyed					C0069415				
3/8 circle reverse cutting needle										
 DS 19 ▼ Easyslide	70 undyed		C0069419	C0069420						
	70 undyed			C0069420N1						
	70 undyed				C0069421					
	70 undyed				C0069421N1					
	70 undyed					C0069422				
1/2 circle reverse cutting needle										
 HS 21 ▼ Easyslide	70 undyed				C0069425					
 HS 23 ▼ Easyslide	70 undyed				C0069430	C0069431	C0069432	C0069433		
 HS 23ss ▼ Easyslide	70 undyed							C0069438		
 HS 26 ▼ Easyslide	70 undyed				C0069440	C0069441				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
 	70 undyed									C0069444
1/2 circle reverse cutting needle										
 	70 undyed						C0069450	C0069451		
 	70 undyed					C0069460				
 	70 undyed							C0069461		
Straight reverse cutting needle										
 	45 undyed					C0069465				
1/2 circle round bodied needle										
 	90 undyed					C0069466	C0069467			
1/2 circle round bodied needle with short cutting point										
 	70 undyed									C0069475
3/8 circle reverse cutting needle										
 	70 undyed				C0069513					
3/8 circle reverse cutting needle with precision point										
 	45 undyed	C0069526	C0069525	C0069524						

NOVOSYN[®]



















ABSORBABLE SUTURES













Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Progressive curved reverse cutting needle with precision point										
 PSMP 16 	45 undyed		C0069602							
1/2 circle reverse cutting needle with precision point										
 HSMP 12 	45 undyed	C0069606								
3/8 circle reverse cutting needle with precision point										
 DSMP 11 	45 undyed	C0069607	C0069608							
 DSMP 13f 	45 undyed	C0069617								
 DSMP 13 	45 undyed	C0069618	C0069619	C0069620						
5/8 circle reverse cutting needle										
 FS 25 	70 undyed				C0069622					
3/8 circle reverse cutting needle with precision point										
 DSMP 19 	45 undyed		C0069628	C0069629						
	 45 undyed			C0069629N1						
	45 undyed				C0069630					
	 45 undyed				C0069630N1					
1/2 circle reverse cutting needle with precision point										
 HSMP 15 	45 undyed			C0069636						
3/8 circle reverse cutting needle with precision point										
 DSMP 24 	45 undyed			C0069655	C0069656	C0069657				
 DSMP 39 	70 undyed				C0069659	C0069660				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with trocar point										
 HRT 37 ⊕ Easyslide	90 undyed					C0069661	C0069662	C0069663	C0069664	
 HRT 40s ⊕ Easyslide	70 undyed						C0069672	C0069671	C0069670	
 HRT 40 ⊕ Easyslide	90 undyed							C0069673		
3/8 circle reverse cutting needle with precision point										
 2xDSMP 7 ▽ Easyslide	45 undyed	C0069709								
 2xDSMP 13 ▽ Easyslide	45 undyed		C0069715	C0069716						
1/2 circle round bodied needle										
 HR 22 ● Easyslide	90 undyed					C0069720				
3/8 circle round bodied needle										
 2xDR 12 ● Easyslide	45 undyed	C0069725								
5/8 circle round bodied needle with trocar point										
 FRT 17 ⊕ Easyslide	70 undyed				C0069755					
1/2 circle round bodied needle with protection point										
 HRB 48 ○ Easyslide	90 undyed								C0069772	

NOVOSYN®












ABSORBABLE SUTURES















Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
  HRC 26	90 undyed					C0069778				
1/2 circle round bodied needle										
  HR 10	45 undyed	C0069802								
1/2 circle round bodied needle with protection point										
  HRB 37s	90 undyed						C0069822			
1/2 circle round bodied needle with short cutting point										
  HRC 37	90 undyed				C0069864	C0069865		C0069867		
  HRC 40	90 70 undyed					C0069920	C0069903			
1/2 circle round bodied needle										
  HR 22s	70 undyed				C0069934					
1/2 circle round bodied needle with short cutting point										
  HRC 40s	70 undyed							C0069958		
  HRC 17	take off 4 x 45 violet					C0088025				
  HRC 26	take off 4 x 45 violet					C0088036				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
 HRC 37 take off	4 x 45 violet						C0088285			
1/2 circle round bodied needle										
 HR 13ss take off	45 violet					C0088320	C0088321	C0088322		
 HR 22 take off	45 violet					C0088325				
 HR 17s take off	45 violet					C0088335				
 HR 17 take off	45 violet				C0088340					
 HR 17 take off	4 x 70 violet				C0088613	C0088614				
 HR 22 take off	4 x 70 violet				C0088629	C0088630				
 HR 26 take off	4 x 70 violet				C0088640	C0088641	C0088642			
 HR 26 take off	4 x 70 violet					C0088647				
 HR 30 take off	4 x 70 violet					C0088660				
 HR 37 take off	4 x 70 violet					C0088669				
 HR 37s take off	4 x 70 violet							C0088696		

NOVOSYN®









ABSORBABLE SUTURES












Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 17 Easyslide take off	4 x 45 violet			C0088713	C0088714					
 HR 30ss Easyslide	70 violet						C0088716			
 HR 13 Easyslide take off	4 x 45 violet		C0088718	C0088719						
 HR 22 Easyslide take off	4 x 45 violet			C0088729	C0088730	C0088731				
 HR 22b Easyblack take off	4 x 45 violet				C0088732					
 HR 26 Easyslide take off	4 x 45 violet			C0088740	C0088741	C0088742				
 HR 30 Easyslide take off	4 x 45 violet					C0088761				
 HR 37 Easyslide take off	4 x 45 violet				C0088769	C0088770				
 HR 26s Easyslide take off	4 x 45 violet					C0088790				
 HR 37s Easyslide take off	4 x 45 violet				C0088794	C0088795				
 HR 26ss Easyslide take off	4 x 45 violet					C0088810				

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 HS 26s 	take off 4 x 45 violet									C0088895
3/8 circle reverse cutting needle										
 DS 24 	take off 4 x 45 violet					C0088953				
1/2 circle round bodied needle										
 HR 37 	take off 4 x 45 undyed									C0089001
 HR 17 	take off 4 x 45 undyed			C0089713						
 HR 26 	take off 4 x 45 undyed					C0089741				
 HR 37s 	take off 4 x 45 undyed									C0089795
1/2 circle reverse cutting needle										
 HS 26s 	take off 4 x 45 undyed					C0089894				

NOVOSYN®























ABSORBABLE SUTURES

















			Thread	USP (metric)							
			length cm colour	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
Ligature reels											
1 thread of 250 cm			1 x 250 violet		G0058715	G0058716	G0058717	G0058719	G0058720		
1 thread of 250 cm			1 x 250 undyed		G0059715	G0059716	G0059717				
angle	Needle length mm diameter microns	Thread length cm colour	USP (metric)								
			8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	
3/8 circle round bodied needle micro needle											
135°	DRm 6  140µ	20 35 violet	G0068145 G0068146								
1/4 circle lancet needle with micro-point											
	VLM 8  Easyslide	70 violet			G0068299						
3/8 circle lancet micro needle											
140°	DLM 6  150µ	15 violet	G0068615								
1/4 circle lancet needle with micro-point											
	2xVLM 6  Easyslide	30 violet			G0068616						
3/8 circle lancet micro needle											
137°	DLM 6s  200µ	30 violet	G0068620	G0068621							
140°	2xDLM 5  150µ	30 violet	G0068699								
140°	2xDLM 6  150µ	45 30 violet	G0068701 G0068702								
137°	2xDLM 6s  200µ	30 45 30 violet	G0068710	G0068711 G0068712							

angle	Needle		Thread length cm colour	USP (metric)							
	length mm	diameter microns		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
3/8 circle lancet micro needle											
137°	2xDLm 6s 	200µ	45 violet	G0068713							
1/2 circle lancet micro needle											
	2xHLM 6 	150µ	30 violet	G0068715							
1/4 circle lancet needle with micro-point											
	VLM 10  <small>Easyslide</small>		45 violet			G0068718					
3/8 circle lancet micro needle											
140°	DLm 6 	150µ	30 violet	G0068724							
1/4 circle lancet needle with micro-point											
	2xVLM 8  <small>Easyslide</small>		20 45 violet			G0068733 G0068734	G0068735	G0068736			
	2xVLM 6  <small>Easyslide</small>		45 violet			G0068737					
	2xVLM 10  <small>Easyslide</small>		45 violet			G0068738	G0068739				
1/2 circle lancet micro needle											
	2xHLM 8 	200µ	30 45 violet	G0068740	G0068741						
1/4 circle lancet needle with micro-point											
	2xVLM 8  <small>Easyslide</small>		20 violet			G0068748					
	VLM 8  <small>Easyslide</small>		45 violet			G0068778					
3/8 circle round bodied needle micro needle											
135°	DRm 8 	140µ	30 violet	G0068796							

NOVOSYN®











ABSORBABLE SUTURES










Needle length mm	Thread length cm colour	USP (metric)								
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	
3/8 circle reverse cutting needle with precision point										
 2xDSMP 6 	45 violet		G0068850							
1/2 circle lancet micro needle										
 2xHLM 8 	45 violet		G0068966	G0068967	G0068968					
1/4 circle lancet needle with micro-point										
 VLM 6 	45 undyed				G0069201					
1/2 circle reverse cutting needle with precision point										
 HSMP 9 	45 undyed					G0069593				
 HSMP 12 	45 undyed				G0069632	G0069633				
 HSMP 15 	45 undyed				G0069635					
1/4 circle lancet needle with micro-point										
 2xVLM 8 	45 undyed		G0069734	G0069735						
1/2 circle lancet micro needle										
 2xHLM 8 	45 undyed		G0069740							
1/2 circle round bodied needle										
 HR 17b  take off	8 x 45 violet							M0068153		
 HR 22b  take off	8 x 45 violet					M0068182	M0068188			
 HR 26b  take off	8 x 45 8 x 70 violet						M0068200 M0068201			

Needle length mm	Thread length cm colour	USP (metric)								
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	
1/2 circle round bodied needle										
	8 x 70 violet									M0068295
										
Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
	8 x 70 violet						M0068305			
										
5/8 circle round bodied needle										
	8 x 90 violet					M0088001				
										
1/2 circle round bodied needle										
	8 x 45 violet		M0088011							
	take off									
1/2 circle round bodied needle with short cutting point										
	8 x 45 violet				M0088037					
	take off									
1/2 circle reverse cutting needle										
	8 x 45 violet				M0088081					
	take off									
	8 x 45 violet			M0088090						
	take off									
1/2 circle round bodied needle										
	8 x 45 violet				M0088101	M0088102	M0088103	M0088104		
	take off									

NOVOSYN[®]







ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle										
 HR 43 Easyslide take off	8 x 45 violet							M0088110	M0088111	
 HR 40s Easyslide take off	8 x 45 violet				M0088115	M0088116	M0088117	M0088118		
 HR 37ss Easyslide take off	8 x 45 violet					M0088131	M0088132			
 HR 37s Easyslide take off	8 x 45 violet			M0088135	M0088136	M0088137	M0088139	M0088140		
 HR 37 Easyslide take off	8 x 45 violet			M0088145	M0088146					
 HR 26ss Easyslide take off	8 x 45 violet				M0088160					
 HR 26s Easyslide take off	8 x 45 violet				M0088164	M0088165	M0088166			
 HR 26 Easyslide take off	8 x 45 violet		M0088170	M0088171	M0088172					
 HR 22 Easyslide take off	8 x 45 violet		M0088177	M0088178	M0088180					
 HR 17 Easyslide take off	8 x 45 violet		M0088185	M0088187						
 HR 13 Easyslide take off	8 x 45 violet	M0088190	M0088191							

Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
3/8 circle reverse cutting needle										
 DS 24 take off	8 x 45 violet				M0088192					
1/2 circle round bodied needle										
 HR 48 take off	8 x 70 violet								M0088204	
 HR 40s take off	8 x 70 violet							M0088207	M0088208	
 HR 37s take off	8 x 70 violet				M0088212	M0088213	M0088214	M0088215		
 HR 37 take off	8 x 70 violet			M0088218	M0088219					
 HR 30ss take off	8 x 70 violet							M0088225		
 HR 26s take off	8 x 70 violet					M0088232	M0088233	M0088234		
 HR 26 take off	8 x 70 violet		M0088237	M0088239						
	8 x 45 undyed				M0089117					
 HR 37s take off	8 x 45 undyed				M0089131	M0089132	M0089133			

NOVOSYN[®]

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle										
	8 x 45 undyed					M0089140				
1/2 circle reverse cutting needle										
	8 x 45 undyed				M0089176					
1/2 circle round bodied needle										
	8 x 70 undyed					M0089504				
1/2 circle round bodied needle with protection point										
	8 x 45 undyed					M0089770	M0089771			
	8 x 45 undyed				M0089819	M0089820	M0089821			
	8 x 45 undyed			M0089828						

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Tel. +34 935 86 62 00 | +34 935 88 10 96 | www.bbraun.es

Information according to Regulation (EU) 2017/745 and Council Directive 93/42/EEC respectively:



Manufacturer

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain

The product trademark "Novosyn" is a registered trademark of B. Braun Melsungen AG or its affiliates.

CE 0123

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form, partial or otherwise, is not permitted.

B08702 0422/PDF/2